

## Product datasheet for **MR202547L4V**

### Tmed3 (NM\_025360) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmed3 (NM_025360) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmed3
Synonyms:	1200002G13Rik; AW546672; P24b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025360
ORF Size:	666 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202547).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_025360.1</a>
RefSeq Size:	1324 bp
RefSeq ORF:	666 bp
Locus ID:	66111
UniProt ID:	<a href="#">Q78IS1</a>
Cytogenetics:	9 E3.1
Gene Summary:	Potential role in vesicular protein trafficking, mainly in the early secretory pathway. Contributes to the coupled localization of TMED2 and TMED10 in the cis-Golgi network (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)