

Product datasheet for **MR202555L4V**

Plaur (BC010309) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Plaur (BC010309) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Plaur
Synonyms:	u-PAR, uPAR, Cd87
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC010309
ORF Size:	666 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202555).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC010309
RefSeq Size:	2076 bp
RefSeq ORF:	668 bp
Locus ID:	18793
Cytogenetics:	7 A3
Gene Summary:	Acts as a receptor for urokinase plasminogen activator. Plays a role in localizing and promoting plasmin formation. Mediates the proteolysis-independent signal transduction activation effects of U-PA.[UniProtKB/Swiss-Prot Function]



[View online »](#)