

Product datasheet for MR202559

Gsta4 (NM_010357) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gsta4 (NM_010357) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gsta4
Synonyms:	GST5.7; mGsta4
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR202559 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGATCGCC

ATGGCAGCCAAACCTAAGCTCTACTACTTTAATGGCAGGGGACGGATGGAGTCGATCCGCTGGCTGCTGG
CTGCGGCTGGAGTGGAGTTTGAGGAAGAATTTCTTGAGACAAGGGAACAGTATGAGAAGATGCAAAAGGA
TGGACACCTGCTTTTCGCCAAGTACCCTTGGTTGAAATCGATGGGATGATGCTGACACAGACCAGGGCC
ATCCTCAGCTACCTCGCTGCCAAGTACAACCTTGTATGGGAAGGACCTGAAGGAGAGATCAGGATTGACA
TGTATGCAGATGGCACCCAGGACCTGATGATGATGATTGCCGTGGCTCCATTTAAAACCCCAAGGAAAA
AGAGGAGAGCTATGATTTGATACTGTCAAGAGCTAAAACCCGTTACTTCCAGTGTTTGAAAAGATTTTA
AAAGACCACGGAGAGGCTTTTCTCGTTGGCAACCAGCTCAGTTGGGCAGACATCCAGCTCCTAGAAGCCA
TTTTGATGGTGGAAAGAACTCAGTGCCCTGTACTGTCCGACTTCCTCTGCTGCAGGCATTTAAGACAAG
AATCAGCAACATTCCTACAATTAAGAAGTTCCTGCAACCCGGAAGTCAGAGGAAGCCTCCTCCAGATGGC
CCCTATGTTGAGGTGGTCAGGACTGTCTGAAGTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR202559 protein sequence
 Red=Cloning site Green=Tags(s)

MAAKPKLYYFNRRGRMESIRWLLAAAGVEFEFEEFLETREQYKMQKDGHLLFGQVPLVEIDGMMLTQTRA
 ILSYLAAYNLYGKDLKERVIRIDMYADGTQDLMMIAVAPFKTPKEEESYDLILSRAKTRYFPVFEKIL
 KDHGEAFLVGNQLSWADIQLLEAILMVEELSAPVLSDFPLLQAFKTRISNIPTIKKFLQPGSQRKPPPDG
 PYVEVVRTVLKF


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_010357

ORF Size: 669 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_010357.3
RefSeq Size:	970 bp
RefSeq ORF:	669 bp
Locus ID:	14860
UniProt ID:	P24472
Cytogenetics:	9 E1
MW:	25.6 kDa
Gene Summary:	Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.[UniProtKB/Swiss-Prot Function]