

Product datasheet for **MR202563**

Krtcap3 (BC016516) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Krtcap3 (BC016516) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Krtcap3
Synonyms:	2010001C09Rik; AW321088; Kcp3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR202563 representing BC016516 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCGCGTGGGACTCGCGCTGATCCTGGTGGGCCACGTGAACCTCCTGGTGGGAGCTGTGTTGCATGGTA
 CTGTTCTGCGGCACGTGGCCAACCCCGAGGCGCGGTCACTCCTGAGTACACCACAGCCAATGTCATTTT
 TGTGGTTCTGGGACTACTGAGCGTTTCAGTGGGACTTGTGGCCCTCTGGCATCCAGGAACCTTCTTCGC
 CCACGACTGCACTGGGCCCTGCTGACGCTGGCCTTGGTGAACCTGCTCTTGTCTGCTGCCTGCTCCATGG
 GCCTCCTCCTGGCTGTATCGCTCACTGTTGCCAATGGTGGCCGCGTCTTATTGCTGATTGCCACCCAGG
 ACTGATGGATCCTTCATACCGCTGGACCAGGGGCTGGACACACTGACTGCTCCTTTGATCCCAAGA
 ATTTATGATACAGCTTTGGCTCTCTGGATCCCTTCTCTGTTTCTGCTGCTGAGCTGAGGCTGCTCTTCTG
 GTTACTGTTGTGTGGCTGCACTCACGCTACGGGAATTGGGCCCTGCAGGAAGGAAGGTTGCAGGAACA
 GCTGCAGGAAGTACAGAGCTTGAACCTCCCGAGTGTAAGGAGGAGAACGTGCAGCTGCTCCATGGC
 CGTCAGGACTTCCAGGCCTTACAGAAGACCTGGTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR202563 representing BC016516
 Red=Cloning site Green=Tags(s)


MRVGLALILVGHVNLVGAHLHGTVLRHVANPRGAVTPEYTTANVISVGSGLLSVSVGLVALLASRNLLR
 PRLHWALLTLALVNLLSAACSMGLLLAVSLTVANGRRRIADCHPGLMDPSIPLDQGPHTDCSFDPTR
 IYDTALALWIPSLFMSAAEALSGYCCVAALTLRGIGPCRKEGLQEQLQELTELELPECKRQENVQLLHG
 RQDFQALQKTWV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: BC016516

ORF Size: 666 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC016516.1](#)

RefSeq Size: 1033 bp


RefSeq ORF: 668 bp

Locus ID: 69815

Cytogenetics: 5 B1

MW: 37.8 kDa

Product images:



Circular map for MR202563