

Product datasheet for **MR202589L1V**

Sigmar1 (NM_011014) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sigmar1 (NM_011014) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sigmar1
Synonyms:	mSig; O; Oprs1; Si; Sig1R; sigma1R
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_011014
ORF Size:	672 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202589).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011014.2
RefSeq Size:	1640 bp
RefSeq ORF:	672 bp
Locus ID:	18391
UniProt ID:	O55242
Cytogenetics:	4 A5



[View online »](#)

Gene Summary:

This gene encodes a transmembrane protein located in the endoplasmic reticulum. The encoded protein is a receptor that binds several endogenous ligands, including N,N-dimethyltryptamine, progesterone and pregnenolone and a variety of of non-opiate compounds. The encoded protein plays a role in regulating the activity of ion channels, acting as a chaperone and protecting cells from oxidative stress. In humans, this receptor has been associated with Alzheimer's and Parkinson's diseases, stroke and numerous disease conditions such as depression, pain and addiction. Alternative splicing results in multiple transcript variants encoding different isoforms.[provided by RefSeq, Nov 2013]