

Product datasheet for **MR202601L4V**

Syng_r2 (NM_009304) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Syng _r 2 (NM_009304) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Syng _r 2
Synonyms:	cellugyrin; Clast2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009304
ORF Size:	672 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202601).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009304.1 , NP_033330.1
RefSeq Size:	1527 bp
RefSeq ORF:	675 bp
Locus ID:	20973
UniProt ID:	O55101
Cytogenetics:	11 E2



[View online »](#)

Gene Summary:

May play a role in regulated exocytosis. In neuronal cells, modulates the localization of synaptophysin/SYP into synaptic-like microvesicles and may therefore play a role in the formation and/or the maturation of this vesicles. May also play a role in GLUT4 storage and transport to the plasma membrane.[UniProtKB/Swiss-Prot Function]