

Product datasheet for **MR202604A1V**

Mouse Apcs (NM_011318) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Apcs (NM_011318) AAV Particle
Tag: Myc-DDK
Symbol: Apcs
Synonyms: Sap
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202604 representing NM_011318
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGACAAGCTACTGCTTTGGATGTTTGTCTTACCAGCCTTCTTTCAGAAGCCTTTTGTGACAGACACC
TCAAGAGGAAAGTATTTGTGTTCCCGAGAGAATCTGAACTGATCATGTGAAGCTGATCCCATCTAGA
GAAACCTCTGCAGAAATTTACTGTGTTTCCGAACCTACAGTGACCTTCCCGCTCTCAGAGTCTTTTC
TCCTACAGTGTCAAGGGCAGAGACAATGAGCTACTAATTTATAAGAAAAAGTTGGAGAATACAGCCTAT
ACATCGGACAATCAAAAGTCACAGTCCGTGGTATGGAAGAATACCTTTCTCCAGTACACCTATGTACCAC
TTGGGAGTCCTCCTCTGGCATTGTTGAATTTTGGGTCAATGGAAGCCTTGGGTAAAAAGTCTCTGCAG
AGGGAATACACTGTGAAAGCCCCACCCAGTATAGTCCTGGGACAGGAGCAGGATAACTACGGAGGAGGGT
TTCAAAGGTCACAGTCTTTGTAGGAGAGTTTTCAGATTTATACATGTGGGACTATGTGCTGACCCACACA
AGACATTCTATTTGTGTACAGAGATTCCCTGTCAATCCTAATATTTGAATTGGCAGGCTCTTAACAT
GAAATAAATGGCTACGTAGTCATCAGGCCCGTGTCTGGGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202604 representing NM_011318
Red=Cloning site Green=Tags(s)

MDKLLLMFVFTSLLSEAFQTDLKRKVFVFPRESETHVKLIPHLEKPLQNFTLCFRTYSDLSRSQSLF
SYSVKGRDNELLYKEKVGESLYIGQSKVTVRGMEEYLSPVHLCTTWESSSGIVEFWNGKPWVKSLQ
REYTVKAPPSIVLGQEQDNYGGGFQRSQSFVGEFSDLYMWDYVLTQDILFVYRDSVPNPNILNWQALNY
EINGYVIRPVWD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_011318
ORF Size:	672 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_011318.2 , NP_035448.2
RefSeq Size:	1012 bp
RefSeq ORF:	675 bp
Locus ID:	20219
UniProt ID:	P12246
Cytogenetics:	1 80.33 cM
MW:	26.7 kDa