

Product datasheet for **MR202608L3V**

Prl2c3 (NM_011118) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Prl2c3 (NM_011118) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Prl2c3
Synonyms:	Ghd1; MRP-2; MRP-3; mrp/plf3; Mrpplf3; PLF-2; PLF-3; Plf2; Prl2c4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011118
ORF Size:	675 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202608).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011118.1 , NP_035248.1
RefSeq Size:	868 bp
RefSeq ORF:	675 bp
Locus ID:	18812
UniProt ID:	P04768
Cytogenetics:	13 A1


[View online »](#)