

Product datasheet for **MR202612L3V**

Zc4h2 (NM_001003916) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Zc4h2 (NM_001003916) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zc4h2
Synonyms:	Gm372; mKIAA1166
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001003916
ORF Size:	675 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202612).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001003916.1
RefSeq Size:	2236 bp
RefSeq ORF:	675 bp
Locus ID:	245522
UniProt ID:	Q68FG0
Cytogenetics:	X C3
Gene Summary:	Plays a role in interneurons differentiation. Involved in neuronal development and in neuromuscular junction formation.[UniProtKB/Swiss-Prot Function]



[View online »](#)