

## Product datasheet for **MR202616L3V**

### **Fkbp3 (NM\_013902) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Fkbp3 (NM_013902) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fkbp3
Synonyms:	25kDa; FKBP-3; FKBP25
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013902
ORF Size:	675 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202616).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013902.2</a>
RefSeq Size:	945 bp
RefSeq ORF:	675 bp
Locus ID:	30795
UniProt ID:	<a href="#">Q62446</a>
Cytogenetics:	12 C1



[View online »](#)

**Gene Summary:**

FK506- and rapamycin-binding proteins (FKBPs) constitute a family of receptors for the two immunosuppressants which inhibit T-cell proliferation by arresting two distinct cytoplasmic signal transmission pathways. PPlases accelerate the folding of proteins (By similarity).  
[UniProtKB/Swiss-Prot Function]