

Product datasheet for **MR202621L3V**

Commd5 (NM_025536) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Commd5 (NM_025536) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Commd5
Synonyms:	2310065H03Rik; AI854466; D15Ertd81e; Hcarg
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025536
ORF Size:	675 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202621).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025536.1 , NP_079812.1
RefSeq Size:	998 bp
RefSeq ORF:	675 bp
Locus ID:	66398
UniProt ID:	Q8R395
Cytogenetics:	15 36.28 cM



[View online »](#)

Gene Summary:

May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes. Negatively regulates cell proliferation. Negatively regulates cell cycle G2/M phase transition probably by transactivating p21/CDKN1A through the p53/TP53-independent signaling pathway. Involved in kidney proximal tubule morphogenesis. Down-regulates activation of NF-kappa-B. [UniProtKB/Swiss-Prot Function]