

Product datasheet for **MR202630L3V**

Snrpb2 (NM_021335) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Snrpb2 (NM_021335) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Snrpb2
Synonyms:	2810052G09Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_021335
ORF Size:	678 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202630).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021335.2 , NP_067310.1
RefSeq Size:	1230 bp
RefSeq ORF:	678 bp
Locus ID:	20639
UniProt ID:	Q9CQI7
Cytogenetics:	2 G1
Gene Summary:	Involved in pre-mRNA splicing as component of the spliceosome. Associated with sn-RNP U2, where it contributes to the binding of stem loop IV of U2 snRNA.[UniProtKB/Swiss-Prot Function]



[View online »](#)