

Product datasheet for **MR202669**

Ms4a4c (NM_029499) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ms4a4c (NM_029499) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ms4a4c
Synonyms:	5830413L19Rik; AI449185; Ms4a9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR202669 representing NM_029499 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAAGGACAGGAACAGACCACCATGGCAGTGGTTCCTGGAGGTGCTCCACCTTCAGAGAATTCTGTTA
 TGAAATCACAAATGTGGAACGAGAACAAGGAGAAATCTTGAAGGGGGAACCCAAAGTCCTTGGGGTTGT
 ACAAGTTATGATTGCTCTCATAACCTCAGCTTCGGAATAATAATTTGGCAAATCTATCTTCTGAACCA
 CTCATTTCTGTGGTCTTAATGGCTCCAATTTGGGGACCAATAATGTTTCATTGTCTCAGGATCCCTGTCAA
 TTGCAGCAGGAGTGAAACCTACAAAAGCCTGATCATCAGCAGTCTAACTCTGAACACTATCACCTCTGT
 GTTGGCTGCAACTGCAAGCATAATGGGTGTAGTCAGTGTGGCTGTGGGTTACAGTTCCTGTTTCGGTAT
 AATTATACAATACCAAGGGTTTGGATATTTTGATGTTAATTTTAAATATGCTAGAATTCTGCATTGCTG
 TGTCCATCTCTGCTTTTGGATGTAAAGCTTCTGTTGTAAGTCCAGCGAGGTTCTTGTAGTGCTACCATC
 AAATCCTGCTGTGACTGTGATGGCACCCCTGTGACGCTTCAACCATGCCACCATCAGAACCAAGGG
 AAAATGTTCCAGAAAATGTATATAAGAACCACTCAGAAGAAATAGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR202669 representing NM_029499
 Red=Cloning site Green=Tags(s)


MQGQEQTMAVVPGGAPPSSENSVMKSQMWNNENKEKFLKGEPKVLGVVQVMIALINLSFGIILANLSSEP
 LISVVLMAPIWGPIMFIVSGSLIAAGVKPTKSLIISSLTLNITISVLAATASIMGVSVAVGSQFPFRY
 NYTITKGLDILMLILNMLEFCIAVSISAFGCKASCNSSEVLVLPSPNAVTVMAPVTLQLPPLPSEHQG
 KNVPENVYKNHSEEIV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9059_d03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_029499

ORF Size: 678 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_029499.3](#), [NP_083775.1](#)

RefSeq Size: 1723 bp


RefSeq ORF: 681 bp

Locus ID: 64380

Cytogenetics: 19 A

MW: 24 kDa

Product images:



Circular map for MR202669