

Product datasheet for **MR202676A1V**

Mouse Yeats4 (NM_026570) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Yeats4 (NM_026570) AAV Particle
Tag: Myc-DDK
Symbol: Yeats4
Synonyms: 4930573H17Rik; B230215M10Rik; GAS41; NuBI-1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202676 ORF sequence
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTTCAAGAGAATGGCGGAATTTGGACCTGACTCCGGCGGGAGAGTGAAGGGGGTTACCATCGTTAAGC
CAATCGTTTATGGCAATGTTGCCAGGTACTTTGGAAAGAAGAGAGAAGAAGACGGGCACACTCACCAGTG
GACTGTGTACGTGAAGCCCTACAGAAACGAGGATATGTCAGCATATGTGAAGAAGATCCAGTTTAAATTA
CAGGAAAGCTACGGCAATCCTCTAAGAGTCGTACCAAGCCTCCATATGAAATCACAGAAACAGGATGGG
GTGAATTTGAAATCATCATCAAGATATTTTCATTGACCCCAACGAGAGACCTGTGACTCTGTACCACTT
ACTGAAGCTCTTCCAGTCTGACACCAATGCCATGCTGGGGAAGAAAACAGTGTTTCAGAGTTCTATGAC
GAAATGATATTTCAAGACCCAACGCAATGATGCAGCAGCTGCTTACGACGTCTCGCCAGCTGACATTGG
GAGCCTATAAGCACGAAACAGAGTTTGCAGAACTTGAAGTAAAAACCAGAGAAAAATTAGAAGCTGCCAA
GAAAAAACAAGCTTTGAGATTGCGGAGCTTAAGGAGAGATTAAGCAAGTCGTGAACTATAAACTGT
TTAAAAAATGAAATCAGGAAGCTTGAAGAAGACGATCAGACGAAAGACATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202676 protein sequence
Red=Cloning site **Green**=Tags(s)

MFKRMAEFGPDSSGRVKGVTIVKPIVYGNVARYFGKKREEDGHTHQWTVYVKPYRNEDMSAYVKKIQFKL
HESYGNPLRVVTKPPYEITETGWGEFEIIKIFFIDPNERPVTLYHLLKLFQSDTNAMLGKKTIVSEFYD
EMIFQDPTAMMQQLTTSRQLTLGAYKHETFAELVKTREKLEAAKKTSFEIAELKERLKASRETINC
LKNEIRKLEEDDQTKDI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_026570
ORF Size:	681 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_026570.1
RefSeq Size:	1418 bp
RefSeq ORF:	684 bp
Locus ID:	64050
UniProt ID:	Q9CR11
Cytogenetics:	10 D2
MW:	26.5 kDa