

Product datasheet for **MR202682A1V**

Mouse Frg2f1 (NM_001004178) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Frg2f1 (NM_001004178) AAV Particle
Tag: Myc-DDK
Symbol: Frg2f1
Synonyms: AA415398; MGC73851; MGC100084
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202682 representing NM_001004178
 Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGATTAGAAACACCCCATCCAGGCCCTTCTGTCAACTGCTGGACCTCCCACTGGACCACAGACCA
 CCTCTGTGACGACATTAGGATCAGAACCAAGTAAAAATGAGCAGAGGTCCAGCCAATCAGACAGCAATTC
 CAGAAAGCACACTGAAGAATCAGCTCAGATGGTGGAAACACCCAGGAAAAGAAGACGCTCTTCCAAGGAC
 TCCAGCCTGGTCAGAAATAGGAAATGAGGGCAGCTTGCATGGGAAAAAGAGACACAAGACTTCAGTAAAGA
 AGACTGAGCCAAAAGCCAATGGTCACCAAGGATGGATTAGCCCATGCCCAGAACCTCAGGATACAGACCC
 ACCACCCATTTCGGAAGAGCCTGGTGAGCTTTCTGCGGGGGAAGTCTGAGGAGATTTATAGAGACACAGTC
 CAGATGCAGGCCCAGAAGCTTGCTCCCTGCTCAGGAGCTCTCTCAGCTCAGACAGCTCTCAG
 AGTCTCTGACTGCCATGGTATACAGCCTTCTACAGCATGGCCAACCAGGAGGCTTTGCCTTCCCTGCTGA
 GGGCTGGCTTGTCCAGCCCCCATGGCTGCCCTCAGGAGCTGTCTGGCAAGGAATCTCAGATCCTTTG
 ATGGAGGGAGGAGAGAAGATAACAGATCCTGCCATCCCATCAGACAAGTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202682 protein sequence
 Red=Cloning site Green=Tags(s)

MDSETPHPGPSCQLLDLPTGPQTTSVRLTGSEPSKNEQRSSQSDSNSRKHTEESAQMVETPRKRRRSSKD
 SSLVRIGNEGSLHGKRRHKTSVKKTEPKANGHQDGLAHCPEPQDTPPPIRKSLVSFLRGKSEEIYRDTV
 QMQAQLGSLLTPEQLSQLRQLESALTAMVYSLLQHGQPGGFAPAEGLVPAPMAAPQELSGKESQSPL
 MEGGEKITDPAIPSDKS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001004178
ORF Size:	681 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001004178.1</u> , <u>NP_001004178.1</u>
RefSeq Size:	1919 bp
RefSeq ORF:	684 bp
Locus ID:	433752
Cytogenetics:	4 D2.1
MW:	24.8 kDa