

Product datasheet for **MR202707L3V**

Hsd17b2 (BC012682) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hsd17b2 (BC012682) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hsd17b2
Synonyms:	AI194836; AI194967; AI255511
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC012682
ORF Size:	1146 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202707).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC012682.1
RefSeq Size:	1280 bp
RefSeq ORF:	684 bp
Locus ID:	15486
Cytogenetics:	8 E1
Gene Summary:	Capable of catalyzing the interconversion of testosterone and androstenedione, as well as estradiol and estrone. Also has 20-alpha-HSD activity. Uses NADH while EDH17B3 uses NADPH (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)