

## Product datasheet for **MR202732L4V**

### **Tomm40 (BC038887) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Tomm40 (BC038887) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tomm40
Synonyms:	Tom40, Mom35
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC038887
ORF Size:	684 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202732).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC038887</a>
RefSeq Size:	1506 bp
RefSeq ORF:	686 bp
Locus ID:	53333
Cytogenetics:	7 A3
Gene Summary:	Channel-forming protein essential for import of protein precursors into mitochondria. [UniProtKB/Swiss-Prot Function]



[View online »](#)