

Product datasheet for **MR202738A1V**

Mouse Zbed3 (NM_028106) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Zbed3 (NM_028106) AAV Particle
Tag: Myc-DDK
Symbol: Zbed3
Synonyms: 0610037K01Rik; 2610005H11Rik; AU018975; AU024588
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202738 representing NM_028106
 Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGCGGCCGGAATTTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAGAGTAAAAAGCCCTGAAGATCACCATGGAAGACAGCCGAGGCTTAATGACCCAGCGGAACAAG
 GCGGCCTCTGTCCCGCGCCCTGGGGCCATCATACTCCGAGGCCTGGGGCTACTTCCATCTGGACCCAGC
 TCAGCCTAGGCACCGGATGATGAGCGCCTGGGCCACCTGCCGCTGTGCGGGCTGCAAGTGGGTGGCCTC
 CCCAACTTCCAGATGTGGACGCGGGCGCTGTGCCAGCACCTGAGTGATGTGCACCTGCCGGAGCTGAAGA
 AGAGCGCTGCTCCGAGCTCGCCGACCACCATGCCCTGCCGCGCCGCCAGCCCAACCATGGCTGCCGA
 GGGCGACTGGGCACGCTGTGGAGCAGATGGGTGAGCTGGCCATGCGGGGCAGCCAGCGGGAGCTGGAG
 CTGGAGCGGCGGAGGCGGCCCTAATGCAGGCGGAACCTGAGCTGGAGCGCAAGCGCCAGGCTCTGAAGC
 AGGAGGCACAATCCGTGGAGCAGGAGCGCCACCAGCTACAGGTAGAGCGTGAAGCACTGAGTAAGTGGAT
 AAAGAAGCAGAGCCCCGGTGCACAGGTCCCGGAACCGCCTCTCCCTGCCTCTGCTCCCAAGGAGGAC
 CCTGACATCCATGACAACAACAGTGACAATGATATGGTCACCAAGGTCTGCTG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202738 representing NM_028106
 Red=Cloning site Green=Tags(s)

MKSKKPLKITMEDSRRLNDPAEQGGLCPAPVGPSYSEAWGYFHLDPAQPRHRMMSAWATCRLCGLQVGGL
 PNFQMWTRALCOHLSDVHLPELKKSAA PSSPTTMCPPPPSPTMAAEGDWARLLEQMGE LAMRGSQRELE
 LERREALMQAELELERKQALKQEAQSVEQERHQLQVEREALSKWIKKQSPGAQVPEPPSPLPLLPKED
 PDIHDNNSDNDMVTKVLL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_028106
ORF Size:	684 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_028106.1
RefSeq Size:	1894 bp
RefSeq ORF:	687 bp
Locus ID:	72114
UniProt ID:	Q9D0L1
Cytogenetics:	13 D1
MW:	25.6 kDa