

Product datasheet for **MR202751L4V**

Pafah1b2 (NM_008775) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pafah1b2 (NM_008775) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pafah1b2
Synonyms:	A1747451; AU021353; mus[b]; Pafahb
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008775
ORF Size:	690 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202751).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008775.3 , NP_032801.2
RefSeq Size:	4241 bp
RefSeq ORF:	690 bp
Locus ID:	18475
UniProt ID:	Q61206
Cytogenetics:	9 A5.2
Gene Summary:	Inactivates PAF by removing the acetyl group at the sn-2 position. This is a catalytic subunit. [UniProtKB/Swiss-Prot Function]



[View online »](#)