

## Product datasheet for MR202752L3

### Scnm1 (NM\_027013) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Scnm1 (NM_027013) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Scnm1
Synonyms:	3110001117Rik; Scnm1-ps
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202752).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	NM_027013
ORF Size:	690 bp



[View online »](#)

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_027013.1</a> , <a href="#">NP_081289.1</a>
<b>RefSeq Size:</b>	888 bp
<b>RefSeq ORF:</b>	561 bp
<b>Locus ID:</b>	69269
<b>UniProt ID:</b>	<a href="#">Q8K136</a>
<b>Cytogenetics:</b>	3 40.74 cM
<b>Gene Summary:</b>	Mutations in the voltage-gated sodium channel gene Scn8a lead to neurological problems in mice. For one particular mutation, Scn8amedJ, mice live to adulthood but have tremors and muscle weakness, among other problems, in all strains except those derived from C57BL6 mice. In these strains, the product of the Scnm1 gene (229 aa) partially overcomes the effects of the Scn8amedJ mutation. However, in C57BL6-derived mice, a one nt change in the penultimate exon creates a premature stop codon, truncating the Scnm1 protein at 186 aa. This truncated protein lacks the ability to overcome the effects of the Scn8amedJ mutation, and these mice suffer paralysis and juvenile death. [provided by RefSeq, Jul 2009]

Product images:



Circular map for MR202752L3