

Product datasheet for **MR202757L3V**

Arl6ip4 (NM_144509) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arl6ip4 (NM_144509) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arl6ip4
Synonyms:	AA408210; AA408365; Aip-4; Srp25
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144509
ORF Size:	690 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202757).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144509.2 , NP_653092.1
RefSeq Size:	1182 bp
RefSeq ORF:	690 bp
Locus ID:	65105
UniProt ID:	Q9JIM93
Cytogenetics:	5 F
Gene Summary:	Involved in modulating alternative pre-mRNA splicing with either 5' distal site activation or preferential use of 3' proximal site.[UniProtKB/Swiss-Prot Function]



[View online »](#)