

Product datasheet for **MR202762L4V**

Tmem126b (NM_026734) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmem126b (NM_026734) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmem126b
Synonyms:	1110001A23Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026734
ORF Size:	693 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202762).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026734.1
RefSeq Size:	890 bp
RefSeq ORF:	693 bp
Locus ID:	68472
UniProt ID:	Q9D1R1
Cytogenetics:	7 E1
Gene Summary:	Chaperone protein involved in the assembly of the mitochondrial NADH:ubiquinone oxidoreductase complex (complex I). Participates in constructing the membrane arm of complex I (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)