

Product datasheet for **MR202764A1V**

Mouse Rnf141 (NM_025999) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Rnf141 (NM_025999) AAV Particle
Tag: Myc-DDK
Symbol: Rnf141
Synonyms: 2610110L04Rik; AA792898; AU022812; ZFP36; ZNF230
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202764 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGACAGCAAATTTTCGGATCAGACGCAGTTGGTTATTAACAAGTTACCAGAAAAGGTAGCAAAGCACG
TCACCTTGTTTCGAGAAAGTGGCTCCTTAACCTACGAAGAGTTTCTTGGGAGAGTAGCTGAACCTCAATGA
TGTAACGTCTAAGGTGGCTGCTGGCCAGGAAAAGCATCTTCTCTTTGAGGTACAACCTGGGTCTGACTCC
TCTGCCTTCTGGAAAGTAGTTGTACGGTGGTCTGTACCAAGATTAACAAAAGCAGTGGCATTGTGGAAG
CATCGCGGATCATGAATTTGTACCAAGTTTATTCAGCTCTATAAAGACATCACAAGCCAAGCTGCTGGAGT
GCTGGCTCAGAGCTCCACCTCTGAAGAGCCTGATGAGAACCCTCCTCGGTACATCTTGCCAGGCTAGT
CTATGGATGGGCAGGTGAAGCAGCTGACTGATGAGGAGGAGTGTGTATCTGCATGGACGGGCGAGCTG
ACCTCATCCTGCCCTGTGCTCACAGCTTCTGTCAGAAGTGCATTGATAAATGGAGTGATCGACACAGGAA
CTGCCCCATTTGTCGCTACAGATGACTGGAGCAAATGAATCATGGGTGGTGTGATGCACCCACTGAA
GATGACATGGCTAACTATATTCTTAACATGGCTGATGAGGCAGGCCAGCCCCACAGGCCA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202764 protein sequence
Red=Cloning site Green=Tags(s)

MGQQISDQTQLVINKLPEKVAKHVTLVRESGLTYEFLGRVAELNDVTAKVAAGQEKHLLFEVQPGSDS
SAFWKVVVRVCTKINKSSGIVEASRIMNLYQFIQLYKDITSQAAGVLAQSSTSEEPDENPSSVTSQAS
LWMGRVKQLTDEEECCICMDGRADLILPCAHSFCQKCIDKWSDRHRNCPICRLQMTGANESWVSDAPTE
DDMANYILNMADEAGQPHRP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_025999
ORF Size:	690 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_025999.2
RefSeq Size:	952 bp
RefSeq ORF:	693 bp
Locus ID:	67150
UniProt ID:	Q99MB7
Cytogenetics:	7 E3
MW:	25.5 kDa