

## Product datasheet for **MR202792L4V**

### Csn2 (NM\_009972) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Csn2 (NM_009972) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Csn2
Synonyms:	Csnb
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009972
ORF Size:	696 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202792).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009972.1</a>
RefSeq Size:	1129 bp
RefSeq ORF:	696 bp
Locus ID:	12991
UniProt ID:	<a href="#">P10598</a>
Cytogenetics:	5 43.56 cM
Gene Summary:	Important role in determination of the surface properties of the casein micelles. [UniProtKB/Swiss-Prot Function]



[View online »](#)