

Product datasheet for **MR202799L4V**

Pafah1b3 (NM_008776) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pafah1b3 (NM_008776) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pafah1b3
Synonyms:	mus[g]; Pafahg
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008776
ORF Size:	699 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202799).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008776.1
RefSeq Size:	892 bp
RefSeq ORF:	699 bp
Locus ID:	18476
UniProt ID:	Q61205
Cytogenetics:	7 A3
Gene Summary:	Inactivates paf by removing the acetyl group at the sn-2 position. This is a catalytic subunit. Plays an important role during the development of brain.[UniProtKB/Swiss-Prot Function]



[View online »](#)