

Product datasheet for **MR202801**

Rassf3 (NM_138956) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rassf3 (NM_138956) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rassf3
Synonyms: AW212023; AW322379
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR202801 representing NM_138956
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCAGCGGCTACAGCAGCCTCGAGGAGGACGAGGACTTTTTCTTCACCGCCAGGACCTCCTTCTTCA
GGAGAGCGCCTCCGGCAAGTCCCGTTCGGGTACCCGGATGTGGAGAAAGAGAAGGAAACGCACAATTA
CCTCAGCAAAGAGGAAATCAAAGAGAAAGTTCATAAATACAATTCAGCCGTTACAGACAAGCTGAAGATG
ACCTTAAATCCAATGGGATCTACTGGCTTCATTAAGGTGCAGATGGAACCTGCAAACCTGCCCAGC
CCTCTCCGAGCCAGCAGCGGTGGCTGCATGAACACGCTTCACATCAGCAGTACAAACACCGTGGGGGA
AGTGATCGAGGCTCTGCTCAGGAAGTTTCTGGTGACCGAGAGTCCACCAAGTTTGCGCTTTATAAGCGT
TGTACCGGGAAGACCAAGTGTATGCCTGCAAGCTCTCCGACCGGAGCACCCGCTGTACCTACGTTTGG
TAGCAGGGCCAGAACCGACACCCTTAGTTTTGTTCTTCGAGAGCATGAAATTGGAGAGTGGGAAGCTTT
CAGCCTACCGGAGCTGCAGAACTTCTGCGCATCTTGACAAGGAAGAGGATGAGCAGCTGCAGAGCCTG
AAGCGGCGCTACACAGCCTACAGGCAGAAGCTGGAGGAAGCCCTGGGCGAGGTGTGGAAGCCCGG

ACGCGTACGCGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR202801 representing NM_138956
Red=Cloning site Green=Tags(s)


MSSGYSSLEEDEDFFFTARTSFFRRAPPGKSRSGQPDVEKEKETHNYLSKEEIKEKVHKYNSAVTDKLM
 TLNSNGIYTGFIKVMELCKPAQSPSESSGGCMNTLHISSTNTVGEVIEALLRKFLVTESPTKFALYKR
 CHREDQVYACKLSDREHPLYLRLVAGPRTDLSFVLRHEHEIGEWAFSLPELQNFLRILDKEEDEQLQSL
 KRRYTAYRQKLEALGEVWPKG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9063_h03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_138956

ORF Size: 696 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_138956.3](#), [NP_620406.1](#)

RefSeq Size: 3425 bp

RefSeq ORF: 699 bp


Locus ID: 192678

UniProt ID: [Q99P51](#)

Cytogenetics: 10 D2

MW: 26.7 kDa

Product images:



Circular map for MR202801