

## Product datasheet for **MR202822L3V**

### **Bcl10 (NM\_009740) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Bcl10 (NM_009740) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bcl10
Synonyms:	A1132454; BCL-10; C81403; CARMEN; cE10; CIPER; CLAP; ME10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009740
ORF Size:	702 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202822).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009740.1</a>
RefSeq Size:	1959 bp
RefSeq ORF:	702 bp
Locus ID:	12042
UniProt ID:	<a href="#">Q9Z0H7</a>
Cytogenetics:	3 H2



[View online »](#)

**Gene Summary:**

Involved in adaptive immune response. Promotes apoptosis, pro-caspase-9 maturation and activation of NF-kappa-B via NIK and IKK. May be an adapter protein between upstream TNFR1-TRADD-RIP complex and the downstream NIK-IKK-IKAP complex. Is a substrate for MALT1.[UniProtKB/Swiss-Prot Function]