

## Product datasheet for **MR202823L3V**

### Tk1 (BC091758) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tk1 (BC091758) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tk1
Synonyms:	Tk1a, Tk1b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC091758
ORF Size:	699 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202823).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC091758</a> , <a href="#">AAH91758</a>
RefSeq Size:	1209 bp
RefSeq ORF:	701 bp
Locus ID:	21877
Cytogenetics:	11 82.96 cM



[View online »](#)