

Product datasheet for **MR202833L4V**

D2Ertd750e (BC031709) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	D2Ertd750e (BC031709) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	D2Ertd750e
Synonyms:	1700025D04Rik; C15orf23; D2Ertd750e; SKAP; Traf4af1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC031709
ORF Size:	702 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202833).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC031709 , AAH31709
RefSeq Size:	1218 bp
RefSeq ORF:	704 bp
Locus ID:	51944
Cytogenetics:	2 E5



[View online »](#)

Gene Summary:

Essential component of the mitotic spindle required for faithful chromosome segregation and progression into anaphase. Promotes the metaphase-to-anaphase transition and is required for chromosome alignment, normal timing of sister chromatid segregation, and maintenance of spindle pole architecture. The astrin (SPAG5)-kinastrin (SKAP) complex promotes stable microtubule-kinetochore attachments. Required for kinetochore oscillations and dynamics of microtubule plus-ends during live cell mitosis, possibly by forming a link between spindle microtubule plus-ends and mitotic chromosomes to achieve faithful cell division.

[UniProtKB/Swiss-Prot Function]