

Product datasheet for **MR202856L3V**

Pomc (NM_008895) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pomc (NM_008895) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pomc
Synonyms:	ACT; ACTH; alp; alph; alpha-MSH; alphaMSH; BE; Beta-LPH; beta-M; beta-MSH; Clip; gamma-Gamma-LPH; gamma-MSH; Npp; PO; Pomc-1; Pomc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008895
ORF Size:	708 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202856).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008895.2
RefSeq Size:	1038 bp
RefSeq ORF:	708 bp
Locus ID:	18976
UniProt ID:	P01193
Cytogenetics:	12 1.99 cM



[View online »](#)

Gene Summary:

This gene encodes a polypeptide hormone precursor that undergoes extensive, tissue-specific, post-translational processing. Processing yields several biologically active peptides, which are involved in diverse cellular functions, such as energy homeostasis, steroidogenesis, and increased melanin production in melanocytes. In mouse deficiency of this gene is associated with obesity, defects in adrenal development, and altered pigmentation. A pseudogene of this gene is located on chromosome 19. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2013]