

Product datasheet for **MR202863L3V**

Trmt11 (BC048703) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Trmt11 (BC048703) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Trmt11
Synonyms:	2410075D05Rik; 3110045I18Rik; AW213713
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC048703
ORF Size:	705 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202863).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC048703 , AAH48703
RefSeq Size:	1574 bp
RefSeq ORF:	707 bp
Locus ID:	73681
Cytogenetics:	10 A4
Gene Summary:	Catalytic subunit of an S-adenosyl-L-methionine-dependent tRNA methyltransferase complex that mediates the methylation of the guanosine nucleotide at position 10 (m2G10) in tRNAs. [UniProtKB/Swiss-Prot Function]



[View online »](#)