

Product datasheet for **MR202873L3V**

Rage (BC087925) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rage (BC087925) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rage
Synonyms:	MGC46883, MOK, RAGE1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC087925
ORF Size:	708 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202873).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC087925.1
RefSeq Size:	1161 bp
RefSeq ORF:	710 bp
Locus ID:	26448
Cytogenetics:	12 F1
Gene Summary:	Able to phosphorylate several exogenous substrates and to undergo autophosphorylation (PubMed:10421840). Negatively regulates cilium length in a cAMP and mTORC1 signaling-dependent manner (PubMed:25243405).[UniProtKB/Swiss-Prot Function]



[View online »](#)