

Product datasheet for **MR202882L4V**

Stx8 (NM_018768) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Stx8 (NM_018768) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Stx8
Synonyms:	0610007H08Rik; 1110002H11Rik; 4930571E13Rik; CARB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018768
ORF Size:	711 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202882).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018768.2
RefSeq Size:	4123 bp
RefSeq ORF:	711 bp
Locus ID:	55943
UniProt ID:	O88983
Cytogenetics:	11 B3
Gene Summary:	Vesicle trafficking protein that functions in the early secretory pathway, possibly by mediating retrograde transport from cis-Golgi membranes to the ER.[UniProtKB/Swiss-Prot Function]



[View online »](#)