

Product datasheet for **MR202893L3V**

Hdgf (NM_008231) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hdgf (NM_008231) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hdgf
Synonyms:	A1118077; D3Ertd299e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008231
ORF Size:	714 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202893).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008231.2
RefSeq Size:	2257 bp
RefSeq ORF:	714 bp
Locus ID:	15191
UniProt ID:	P51859
Cytogenetics:	3 38.78 cM
Gene Summary:	Heparin-binding protein, with mitogenic activity for fibroblasts. Acts as a transcriptional repressor (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)