

## Product datasheet for **MR202895L4V**

### **Cdca4 (NM\_028023) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Cdca4 (NM_028023) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cdca4
Synonyms:	2410018C03Rik; HEPP; SEI-3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_028023
ORF Size:	714 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202895).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_028023.3</a> , <a href="#">NP_082299.1</a>
RefSeq Size:	2432 bp
RefSeq ORF:	714 bp
Locus ID:	71963
UniProt ID:	<a href="#">Q9CWM2</a>
Cytogenetics:	12 F1
Gene Summary:	May participate in the regulation of cell proliferation through the E2F/RB pathway (By similarity). May be involved in molecular regulation of hematopoietic stem cells and progenitor cell lineage commitment and differentiation.[UniProtKB/Swiss-Prot Function]



[View online »](#)