

## Product datasheet for **MR202930L4V**

### Ube2r2 (NM\_026275) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ube2r2 (NM_026275) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ube2r2
Synonyms:	1200003M11Rik; Cdc34b; Ubc3b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026275
ORF Size:	717 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202930).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026275.4</a> , <a href="#">NP_080551.1</a>
RefSeq Size:	3312 bp
RefSeq ORF:	717 bp
Locus ID:	67615
UniProt ID:	<a href="#">Q6ZWZ2</a>
Cytogenetics:	4 A5
Gene Summary:	Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro catalyzes monoubiquitination and 'Lys-48'-linked polyubiquitination. May be involved in degradation of katenin.[UniProtKB/Swiss-Prot Function]



[View online »](#)