

Product datasheet for MR202937

Psme2 (NM_011190) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Psme2 (NM_011190) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Psme2
Synonyms:	AA589371; AI788882; PA28b
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR202937 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCAAGCCTTGTTGGGTCGCTTGAGTGGAGAAGCCCGAAAACAGGTGGATGTCTTCAGACAAAATC
TTTTCCAGGAGGCTGATGACTTCCTCTGCACTTTCTTGCCACGGAAAATCATATCCCTGAGTCAGCTCTT
GCAGGAGGATTCCTCAATGTGGCCGACCTCTCCTCCCTCCGGGCTCCTCTGGACATCCCTATCCCAGAT
CCTCCACCCAAGGATGACGAGATGGAAACAGACAAGCAGGAGAAGAAGGAAGTCCCTAAGTGTGGCTACC
TCCCGGGGAACGAGAAGCTTCTGGCCTTGCTCGCTTTGGTTAAGCCAGAAGTCTGGACTCTCAAAGAGAA
GTGCAATTCGTAATCACATGGATCCAGCACCTGATCCCCAAGATTGAGGATGGAAATGATTTTGGGGTG
GCAATTCAGGAGAAGGTGCTGGAGAGAGTGAATGCTGTCAAGACCAAAGTAGAAGCTTTCAGACAACCA
TTTCCAAGTACTTCTCAGAACGAGGGGACGCTGTAGCCAAGGCTTCCAAGGACACCCATGTAATGGATTA
CCGGGCCCTTGGTGCATGAGCGAGATGAAGCAGCCTATGGGGCGCTTAGGGCCATGGTGTGGACCTGAGG
GCCTTCTATGCTGAGCTTTATCATATCATCAGCAGCAACCTAGAGAAAATCGTCAACCCGAAGGGTGAAG
AAAAGCCATCGATGTAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR202937 protein sequence
 Red=Cloning site Green=Tags(s)

MAKPCGVRLSGEARKQVDVFRQNLFQEADDFLCTFLPRKIISLSQLLQEDSLNVADLSSLRAPLDIIPD
 PPPKDDMETDKQEKKEVPKCGYLPGNEKLLALLALVKPEVWTLKEKILVITWQHLPKIEDGNDFGV
 AIQEKVLERNVAVKTKVEAFQTTISKYFSEGRDAVAKASKDTHVMDYRALVHERDEAAYGALRAMVLDLR
 AFYAELYHISSNLEKIVNPKGEEKPSMY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_011190

ORF Size: 717 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_011190.3](#), [NP_035320.1](#)

RefSeq Size: 906 bp

RefSeq ORF: 720 bp


Locus ID: 19188

UniProt ID: [P97372](#)

Cytogenetics: 14 28.19 cM

MW: 27.1 kDa

Product images:



Circular map for MR202937