

Product datasheet for **MR202939A1V**

Mouse Rnf151 (NM_026205) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Rnf151 (NM_026205) AAV Particle
Tag:	Myc-DDK
Symbol:	Rnf151
Synonyms:	1700010O16Rik; AV047371
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR202939 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGTGGTGGGTACGATCTCAACCTCTTTGCCAGCCCTCCTGACTGCAAGTTCCTGTGTTCTGTCTGCC
ATGGGGTTCTCAAGAGGCCAACGAGGTTGCCGTGCAGTCACATCTTCTGTAAAAAGTGCATCTTCAGGTG
GCTAGCCAGACAAAACACCTGTCCATGCTGTAGAAAAGAGGTGACAAGGAGAAAGATGGTCGAAGTGAAT
AAACTCCGAAAAACCATTGGCCGCCTACAAGTCAAGTGCAAGAACGCTGCAGCCGGCTGCTTGGACACAC
ACCCCTCTGGCTCATCGCAAGGAGCATCAGGACTCATGCCCTTTGAGCTGATGGCCTGCCCAATGAGGG
TTGCACTGTGCAGGTGCTGCGAGGGGTCTGGATGAACACCGCCAGCATTGCCAGCAGAATGGACAGCAG
CGCTGCCCTTAGGCTGTGGGTCCACCCTGGCTGCCTTGAAGGTGAGCACCACAACCTGTACCGCGAGC
TACGTGATGCTTGGGTCCAACGCCATGAGCGCAACCGGACCCTGCTGCTGGGCCTGCTGGGGCGTGACG
CCGGGTGCACCTCACCACCAGCATCATCCACCAGCAGCTGGCAGAGCTGAGCAACTTCCTGGAGGATGAT
GATAACCTGCTCCTGAATGCCAGGTGCAAGAGACTGAGGTTACCCCGGAGGCTGAGATGAGGGGTACGC
AGGGCCAAAGTGCTTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR202939 protein sequence Red=Cloning site Green=Tags(s) MSGGYDLNLFASPPDCKFLCSVCHGVLKRPTLPCSHIFCKKCIFRWLARQNTCPCCRKEVTRRKMEVN KLRKTIGRLQVKCKNAAAGCLDTHPLAHRKEHQDSCPFELMACPNEGCTVQVLRGVLDEHRQHCQNGQQ RCPLGCGSTLAALERGEHHNCYRELDAWVQRHERNRTLLLGLLGRVRRVHLTTSIIHQQLAQLSNFLEDD DNLLLNQVQETEVTPAEEMRGTOGQSVL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_026205
ORF Size:	717 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_026205.3, NP_080481.1</u>
RefSeq Size:	1189 bp
RefSeq ORF:	720 bp
Locus ID:	67504
UniProt ID:	<u>Q9CQ29</u>
Cytogenetics:	17 A3.3
MW:	27.2 kDa