

Product datasheet for **MR202946**

Egln3 (NM_028133) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Egln3 (NM_028133) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Egln3
Synonyms:	2610021G09Rik; AI505553; AI648162; Hif-p4h-3; Phd3; SM-20
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR202946 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCTCTGGGACACATCATGAGGCTGGATCTGGAGAAGATCGCCCTGGAGTACATCGTGCCCTGTCTGC
ACGAGGTCGGCTTCTGCTACCTGGACAACCTCCTAGGCGAGGTGGTGGGCGACTGCGTGCTGGAGCGAGT
CAAGCAGCTGCACTACAACGGGGCCCTGCGTGACGGCCAGCTGGCCGGCCGCGCTGGCGTCTCCAAG
CGGCACCTGCGAGGCGACCAGATCACGTGGATCGGGGCAACGAGGAGGGCTGCGAGGCCATCAACTTCC
TCCTGTCCCTCATCGACCGGCTGGTCTGTACTGCGGGAGCCGGCTGGGCAAATACTATGTCAAGGAGCG
GTCCAAGGCAATGGTGGCTTGTATCCAGGAAATGGGACAGGTTATGTTGCCATGTGGACAACCCCAAT
GGTGATGGCCGCTGTATCACCTGTATCTACTACCTGAATAAGAATTGGGACGCCAAGTTACACGGAGGGG
TCCTGCGGATATTTCCGGAAGGGAAATCGTTTGTAGCAGACGTGGAGCCATTTTGGACAGACTTCTGTT
CTTCTGGTCAGACCGCAGGAATCCACATGAAGTCCAGCCCTCCTATGCCACCAGGTACGCTATGACTGTC
TGGTACTTCGATGCTGAAGAAAGGGCAGAAGCCAAAAAGAAATTCAGGAATTTAACTAGAAAACTGAAT
CTGCTCTTGCTAAAGAC

ACGCGTACGCGGCGGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR202946 protein sequence
 Red=Cloning site Green=Tags(s)

MPLGHIMRLDLEKIALEYIVPCLHEVGFVFCYLDNFLGEVVGDCVLERVKQLHYNGALRDGQLAGPRAGVSK
 RHLRGDQITWIGNEEGCEAINFLLSLIDRLVLYCGSRLGKYYVKERSKAMVACYPGNGTGYVRHVDNPN
 GDGRCITCIYYLNKNWDAKLHGGVLRIFPEGKSFVADVEPIFDRLFFWSDRRNPHEVQPSYATRYAMTV
 WYFDAEERAEAKKKFRNLTRKTESALAKD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_028133

ORF Size: 720 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_028133.1](#), [NP_082409.1](#)

RefSeq Size: 2693 bp

RefSeq ORF: 720 bp

Locus ID: 112407


UniProt ID: [Q91UZ4](#)

Cytogenetics: 12 C1

MW: 27.3 kDa

Gene Summary: Plays a crucial role in DNA damage response (DDR) by hydroxylating TELO2, promoting its interaction with ATR which is required for activation of the ATR/CHK1/p53 pathway (By similarity). Cellular oxygen sensor that catalyzes, under normoxic conditions, the post-translational formation of 4-hydroxyproline in hypoxia-inducible factor (HIF) alpha proteins. Hydroxylates a specific proline found in each of the oxygen-dependent degradation (ODD) domains (N-terminal, NODD, and C-terminal, CODD) of HIF1A. Also hydroxylates HIF2A. Has a preference for the CODD site for both HIF1A and HIF2A. Hydroxylation on the NODD site by EGLN3 appears to require prior hydroxylation on the CODD site. Hydroxylated HIFs are then targeted for proteasomal degradation via the von Hippel-Lindau ubiquitination complex. Under hypoxic conditions, the hydroxylation reaction is attenuated allowing HIFs to escape degradation resulting in their translocation to the nucleus, heterodimerization with HIF1B, and increased expression of hypoxia-inducible genes. ELGN3 is the most important isozyme in limiting physiological activation of HIFs (particularly HIF2A) in hypoxia. Also hydroxylates PKM in hypoxia, limiting glycolysis. Under normoxia, hydroxylates and regulates the stability of ADRB2. Regulator of cardiomyocyte and neuronal apoptosis. In cardiomyocytes, inhibits the anti-apoptotic effect of BCL2 by disrupting the BAX-BCL2 complex. In neurons, has a NGF-induced proapoptotic effect, probably through regulating CASP3 activity. Also essential for hypoxic regulation of neutrophilic inflammation. Target proteins are preferentially recognized via a LXXLAP motif.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR202946