

Product datasheet for **MR202963L3V**

Tpmt (NM_016785) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tpmt (NM_016785) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tpmt
Synonyms:	AW106912
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016785
ORF Size:	723 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202963).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016785.1 , NP_058065.1
RefSeq Size:	2715 bp
RefSeq ORF:	723 bp
Locus ID:	22017
UniProt ID:	O55060
Cytogenetics:	13 24.5 cM



[View online »](#)

Gene Summary:

Catalyzes the S-methylation of thiopurine drugs such as 6-mercaptopurine (also called mercaptopurine, 6-MP or its brand name Purinethol) using S-adenosyl-L-methionine as the methyl donor (PubMed:18484748). TPMT activity modulates the cytotoxic effects of thiopurine prodrugs. A natural substrate for this enzyme has yet to be identified. [UniProtKB/Swiss-Prot Function]