

Product datasheet for **MR202966L4V**

Efhd1 (NM_028889) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Efhd1 (NM_028889) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Efhd1
Synonyms:	4931430I01Rik; A1452351; mitocalcin; PP3051
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_028889
ORF Size:	723 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202966).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_028889.2 , NP_083165.1
RefSeq Size:	1805 bp
RefSeq ORF:	723 bp
Locus ID:	98363
UniProt ID:	Q9D4J1
Cytogenetics:	1 C5
Gene Summary:	Acts as a calcium sensor for mitochondrial flash (mitoflash) activation, an event characterized by stochastic bursts of superoxide production (By similarity). May play a role in neuronal differentiation (PubMed:16336229).[UniProtKB/Swiss-Prot Function]



[View online »](#)