

Product datasheet for **MR202977L4V**

Gstt3 (NM_133994) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gstt3 (NM_133994) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gstt3
Synonyms:	A1118089
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133994
ORF Size:	726 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202977).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133994.3 , NP_598755.1
RefSeq Size:	1873 bp
RefSeq ORF:	726 bp
Locus ID:	103140
UniProt ID:	Q99L20
Cytogenetics:	10 C1



[View online »](#)

Gene Summary:

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Shows high activity towards 4-nitrobenzyl chloride (4-NBC). Also has lower activity towards 1,2-epoxy-3-(p-nitrophenoxy)propane (EPNP), cumene hydroperoxide, 1-chloro-2,4-dinitrobenzene (CDNB), 7-chloro-4-nitrobenzo-2-oxa-1,3-diazole (NBD-Cl), and ethacrynic acid.[UniProtKB/Swiss-Prot Function]