

Product datasheet for **MR202986L4V**

Qdpr (NM_024236) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Qdpr (NM_024236) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Qdpr
Synonyms:	2610008L04Rik; D5Ertd371e; Dhpr; PKU2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024236
ORF Size:	726 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202986).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024236.2 , NP_077198.1
RefSeq Size:	1394 bp
RefSeq ORF:	726 bp
Locus ID:	110391
UniProt ID:	Q8BVI4
Cytogenetics:	5 24.9 cM
Gene Summary:	The product of this enzyme, tetrahydrobiopterin (BH-4), is an essential cofactor for phenylalanine, tyrosine, and tryptophan hydroxylases.[UniProtKB/Swiss-Prot Function]



[View online »](#)