

Product datasheet for **MR202987L3V**

Tnfsf13 (NM_001159505) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tnfsf13 (NM_001159505) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tnfsf13
Synonyms:	2310026N09Rik; April; Tall2; Tnlg7b; Trdl1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001159505
ORF Size:	726 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202987).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001159505.1 , NP_001152977.1
RefSeq Size:	1676 bp
RefSeq ORF:	723 bp
Locus ID:	69583
UniProt ID:	Q9D777
Cytogenetics:	11 B3
Gene Summary:	Cytokine that binds to TNFRSF13B/TACI and to TNFRSF17/BCMA. Plays a role in the regulation of tumor cell growth. May be involved in monocyte/macrophage-mediated immunological processes.[UniProtKB/Swiss-Prot Function]



[View online »](#)