

## Product datasheet for MR203006L4

### Pbx1 (BC002244) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Pbx1 (BC002244) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Pbx1
Synonyms:	2310056B04Rik; D230003C07Rik; Pbx-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203006).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	BC002244
ORF Size:	726 bp



[View online »](#)

<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	<p>This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.</p>
<b>Components:</b>	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<p><a href="#">BC002244</a>, <a href="#">AAH02244</a></p>
<b>RefSeq Size:</b>	<p>1491 bp</p>
<b>RefSeq ORF:</b>	<p>728 bp</p>
<b>Locus ID:</b>	<p>18514</p>
<b>Cytogenetics:</b>	<p>1 75.95 cM</p>
<b>Gene Summary:</b>	<p>This gene encodes a homeobox protein that belongs to the three-amino-acid loop extension/Pre-B cell leukemia transcription factor (TALE/PBX) family of proteins. The encoded protein is involved in several biological processes during embryogenesis including steroidogenesis, sexual development and the maintenance of hematopoietic stem cells. This protein functions in the development of several organ systems and plays a role in skeletal patterning and programming. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]</p>

Product images:



Circular map for MR203006L4