

Product datasheet for **MR203026L3V**

LOC629446 (BC057932) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Symbol: LOC629446

Mammalian Cell Selection: Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK

ACCN: BC057932

ORF Size: 729 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR203026).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [BC057932](#)

RefSeq Size: 1558 bp

RefSeq ORF: 731 bp

Locus ID: 629446

Gene Summary: DISCONTINUED: This record has been withdrawn by NCBI. This record was defined by BC057932.1 which represents retroviral sequence.

