

Product datasheet for **MR203050L4V**

C14orf166 (Rtraf) (NM_026528) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C14orf166 (Rtraf) (NM_026528) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rtraf
Synonyms:	2700060E02Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026528
ORF Size:	735 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203050).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026528.3 , NP_080804.1
RefSeq Size:	1510 bp
RefSeq ORF:	735 bp
Locus ID:	68045
UniProt ID:	Q9CQE8
Cytogenetics:	14 A3
Gene Summary:	RNA-binding protein involved in modulation of mRNA transcription by Polymerase II. Component of the tRNA-splicing ligase complex and is required for tRNA ligation. May be required for RNA transport.[UniProtKB/Swiss-Prot Function]



[View online »](#)