

Product datasheet for **MR203060L3V**

Cnpy4 (NM_178612) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cnpy4 (NM_178612) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cnpy4
Synonyms:	2610019P18Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_178612
ORF Size:	738 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203060).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178612.3 , NP_848727.1
RefSeq Size:	1781 bp
RefSeq ORF:	738 bp
Locus ID:	66455
UniProt ID:	Q8BQ47
Cytogenetics:	5 G2
Gene Summary:	Plays a role in the regulation of the cell surface expression of TLR4.[UniProtKB/Swiss-Prot Function]



[View online »](#)