

Product datasheet for **MR203061A1V**

Mouse Fgfbp3 (NM_028263) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Fgfbp3 (NM_028263) AAV Particle
Tag:	Myc-DDK
Symbol:	Fgfbp3
Synonyms:	2610306H15Rik; Fgf-bp3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR203061 representing NM_028263 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGTCCTCCGAGGCCGCGGGCGTCTCTCCCGCTGACGCTGCTGTTGCTGCTGGGAGGCTGCCTCC
TCTCGGCTGCCGGGAGGGACAAGGGGGCAGCAGGCAGGGAAGTGACCCGGGCCTCCCGACCCACGGTCGG
CTCTCCGGTCGCTTCGTGAGTCCAGAGCAGCACGCTGCAGCTGGCAGCTCTTGGTACCCGCCCGGGG
ACGCCGACTGGCGGAGAGCTGGCCTTGCCTGCCAGACCCGGCGGGCGGAGTCTGCATTGCGCCTACC
GCGGGCATCCGGAGCGCTGCGCGGCCACCGCGCCCGCGCGCGCACTACTGGAGACGGCTGCTGGGCGC
GCTGCGCAGGAGACCACGGCCCTGCCTGGACCCCGCGCCGCTGCCGCCACGCCTGTGCGCCCGCAAGACG
GCTGGCTCAGACCTGCACTCGCCCGCTACCCAGCCTACCCGCGCGTCCAGCGAGCCGCCCGGTCCA
GGGCGCGGAGCCCGGCTCGGTGCGCAGTGGTGGGTCCCCAAGTCCAGCCGAAAGAAGCCCTT
GCTAGTGAAGTCCAACCTCGGCGGGAGGAAAGCGGGCTCAGACCCAGTTCCGAGCCTCCTGCGGCGGCC
GGATTCCAACCAATGGGCTGGACCAGAATGCCGAGCTCACTGAGACCTACTGCACCGAGAAGTGGCATT
CCCTTTGCAACTTCTTTGTCAACTTCTGGAATGGC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR203061 representing NM_028263 Red=Cloning site Green=Tags(s) MSPPRPRASLSPLTLLLLLGGCLLSAAGRDKGAAGREVTRASRPTVGSSGRFVSPEQHACSWQLLPAPG TPTGGELALRCQTPGGASLHCAYRGHPERCAATGARRAHYWRRLLGALRRRPRPCLDPAPLPPRLCARKT AGSDLHSPAHPSLPARPSEPSPRSRARSQSVSPSSQPEKKPLLKSNSSGGRKAGSDPVPEPPAAA GFQPNGLDQNAELTETCYCTEKWHSLCNFFVNFVWG SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_028263
ORF Size:	735 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_028263.1</u>
RefSeq Size:	1775 bp
RefSeq ORF:	753 bp
Locus ID:	72514
UniProt ID:	<u>Q1HCM0</u>
Cytogenetics:	19 C2
MW:	27.1 kDa