

## Product datasheet for **MR203075L3V**

### **Bdh2 (NM\_027208) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Bdh2 (NM_027208) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bdh2
Synonyms:	1810026B04Rik; Dhrrs6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_027208
ORF Size:	738 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203075).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_027208.1</a>
RefSeq Size:	1099 bp
RefSeq ORF:	738 bp
Locus ID:	69772
UniProt ID:	<a href="#">Q8JZV9</a>
Cytogenetics:	3 G3



[View online »](#)

**Gene Summary:**

Dehydrogenase that mediates the formation of 2,5-dihydroxybenzoic acid (2,5-DHBA), a siderophore that shares structural similarities with bacterial enterobactin and associates with LCN2, thereby playing a key role in iron assimilation and homeostasis (PubMed:20550936, PubMed:24863067). Plays a role in susceptibility to bacterial infection by providing an assimilable source of iron that is exploited by pathogenic bacteria (PubMed:24863067). Also acts as a 3-hydroxybutyrate dehydrogenase (By similarity).[UniProtKB/Swiss-Prot Function]